

2

AD

TECHNICAL REPORT 8903

DETERMINATION OF DIMETHYL SULFOXIDE (DMSO), ETHANOL (ETOH),
FORMAMIDE (F) AND GLYCEROL/FORMAL (GF) BY HIGH PERFORMANCE LIQUID
CHROMATOGRAPHY (HPLC)

ERNST E. BRUEGGEMANN
ALAN B. ROSENCRANCE

AD-A231 385

DTIC
ELECTE
JAN 23 1991
S B D

U S ARMY BIOMEDICAL RESEARCH & DEVELOPMENT LABORATORY

Fort Detrick

Frederick, MD 21701-5010

January 1989

Approved for public release;
distribution unlimited.



U S ARMY MEDICAL RESEARCH & DEVELOPMENT COMMAND
Fort Detrick
Frederick, MD 21701-5012

REPORT DOCUMENTATION PAGE

Form Approved
OMB No 0704-0188

1a REPORT SECURITY CLASSIFICATION Unclassified			1b RESTRICTIVE MARKINGS		
2a SECURITY CLASSIFICATION AUTHORITY			3 DISTRIBUTION/AVAILABILITY OF REPORT Approved for public release; distribution unlimited.		
2b DECLASSIFICATION/DOWNGRADING SCHEDULE					
4. PERFORMING ORGANIZATION REPORT NUMBER(S) Technical Report #8903			5. MONITORING ORGANIZATION REPORT NUMBER(S)		
6a. NAME OF PERFORMING ORGANIZATION U.S. Army Biomedical Research & Development Laboratory		6b. OFFICE SYMBOL (If applicable) SGRD-UBG-R		7a. NAME OF MONITORING ORGANIZATION	
6c. ADDRESS (City, State, and ZIP Code) Building 568, Fort Detrick Frederick, MD 21702-5010		7b. ADDRESS (City, State, and ZIP Code)			
8a. NAME OF FUNDING/SPONSORING ORGANIZATION U.S. Army Biomed- ical Research & Development Lab		8b. OFFICE SYMBOL (If applicable) SGRD-UBG-R		9 PROCUREMENT INSTRUMENT IDENTIFICATION NUMBER	
8c. ADDRESS (City, State, and ZIP Code) Building 568, Fort Detrick Frederick, MD 21702-5010		10 SOURCE OF FUNDING NUMBERS			
		PROGRAM ELEMENT NO		PROJECT NO	TASK NO.
				WORK UNIT ACCESSION NO	
11. TITLE (Include Security Classification) (U) Determination of Dimethyl Sulfoxide (DMSO), Ethanol, (ETOH), Formamide (F), and Glycerol/ Formal (GF) by High Performance Liquid Chromatography (HPLC).					
12. PERSONAL AUTHOR(S) ERNST E. BRUEGGEMANN and ALAN B. ROSENCRANCE					
13a. TYPE OF REPORT Technical		13b. TIME COVERED FROM 1/88 TO 5/88		14. DATE OF REPORT (Year, Month, Day) 1989 January 30	
				15. PAGE COUNT 13	
16. SUPPLEMENTARY NOTATION					
17. COSATI CODES			18. SUBJECT TERMS (Continue on reverse if necessary and identify by block number)		
FIELD	GROUP	SUB-GROUP	High Performance Liquid Chromatograph (HPLC), Dimethyl/ Sulfoxide (DMSO), Ethanol (ETOH), Formamide (F), Glycerol/ Formal (GF), Analytical Method -		
19. ABSTRACT (Continue on reverse if necessary and identify by block number) Liquid chromatographic (LC) methods have been developed for the determination of dimethyl sulfoxide (DMSO), ethanol (ETOH), formamide (F), and glycerol/formal (GF), in aqueous samples. The samples are directly injected onto a high performance liquid chromatograph (HPLC) without further sample preparation. Separation of DMSO and F was achieved by using a C18 column and a mobile phase containing water/methanol. Separation of GS and ETOH was achieved by using a Dextro-Pak column and a mobile phase consisting of water. A programmable ultraviolet (UV) spectrophotometric detector was used to monitor the LC effluent for DMSO and F. A refractive index detector was used to monitor the LC effluent for ETOH and GF. The precision and accuracy of all methods are included in the report.					
20. DISTRIBUTION/AVAILABILITY OF ABSTRACT <input checked="" type="checkbox"/> UNCLASSIFIED/UNLIMITED <input type="checkbox"/> SAME AS RPT <input type="checkbox"/> DTIC USERS			21 ABSTRACT SECURITY CLASSIFICATION Unclassified		
22a. NAME OF RESPONSIBLE INDIVIDUAL Ernst E. Brueggemann			22b TELEPHONE (Include Area Code) (301) 663-2340		22c. OFFICE SYMBOL SGRD-UBG-R

NOTICE

Disclaimers

The views, opinions, and/or findings contained in this report are those of the author(s) and should not be construed as official Department of the Army position, policy, or decision, unless so designated by other official documentation.

Research was conducted in compliance with the Animal Welfare Act, and other Federal statutes and regulations relating to animals and experiments involving animals and adheres to principles stated in the Guide for the Care and Use of Laboratory Animals, NIH publication 86-23, 1985 edition.

Disposition

Destroy this report when it is no longer needed. Do not return it to the originator.

Accession For	
NTIS GRA&I	<input checked="checked" type="checkbox"/>
DTIC TAB	<input type="checkbox"/>
Unannounced	<input type="checkbox"/>
Justification	
By _____	
Distribution/	
Availability Codes	
Dist	Avail and/or Special
A-1	

TABLE OF CONTENTS

	<u>Page</u>
INTRODUCTION AND OBJECTIVES.....	1
METHODS AND MATERIALS.....	1
Analytical Instrumentation.....	1
Chemicals.....	1
Chemical and Physical Properties.....	2
HPLC Conditions.....	3
Preparation of Stock and Standard Solutions.....	4
Sample Preparation and Handling.....	5
RESULTS.....	5
DISCUSSION.....	5
CONCLUSION.....	6
LITERATURE CITED.....	7
DISTRIBUTION LIST.....	8

TABLES

1. Standard Curve Data.....	9
2. Precision Data for Dimethyl Sulfoxide and Ethanol.....	10
3. Precision Data for Formamide and Glycerol/Formal.....	11
4. Accuracy Data for Dimethyl Sulfoxide and Ethanol.....	12
5. Accuracy Data for Formamide and Glycerol/Formal.....	13

FIGURES

1. HPLC Chromatogram of Dimethyl Sulfoxide in Water.....	2
2. HPLC Chromatogram of Ethanol in Water.....	2
3. HPLC Chromatogram of Formamide in Water.....	2
4. HPLC Chromatogram of Glycerol/Formal in Water.....	3

INTRODUCTION AND OBJECTIVES

Currently, investigations are being performed at the U.S. Army Biomedical Research and Development Laboratory, Research Methods Branch, to determine the teratogenicity of ethanol, dimethyl sulfoxide, formamide and glycerol/formal on Xenopus laevis embryos.

In support of these investigations, it was our task to develop a separate High Performance Liquid Chromatographic (HPLC) method for each of these four compounds in aqueous solutions.

Although other analytical methods exist,¹⁻³ HPLC (employing reverse phase columns) appears to be the analytical method of choice because the water miscible samples can be directly injected onto the HPLC without any prior sample clean-up. In addition, different types of detectors can be easily interchanged in the HPLC system to detect the appropriate compound in solution.

Each HPLC method reported herein is intended to be quick, sensitive and reproducible for the determination of each compound in aqueous solution.

METHODS AND MATERIALS

ANALYTICAL INSTRUMENTATION

A Waters liquid chromatographic system, (Millipore Corp., Waters Chromatography Division, Milford, MA), was employed throughout the study. The system consisted of the following components: a model M6000A solvent delivery system, a model 710B WISP autosampler, a model 721 programmable system controller, a model 730 data module, a model R-400 refractive index detector, and a model 480 LC Spectrophotometer.

CHEMICALS

The methanol used was "HPLC grade" from Burdick and Jackson Laboratories (Muskegon, MI). Reagent grade water was obtained with a Milli-Q System, (Millipore, Bedford, MA), and had a resistance of 18 megohms-cm.


Dimethyl Sulfoxide (99 + % pure), glycerol/formal (98% pure) and formamide (99% pure) were obtained from the Aldrich Chemical Company, Inc., (Milwaukee, WI) and were used without further purification.

The ethanol (95 + % pure) used throughout this study was from U.S. Industrial Chemicals Co. (Tuscola, IL).

CHEMICAL AND PHYSICAL PROPERTIES

DMSO


CAS Registry Number: 67-68-5
RTECS Reference Number: PV6210000
M.W. 78.13
Chemical Formula: C₂H₆OS
Structure:



Synonyms: DEMESO; DEMASORB; DEMAVET; DERMASORB; DIMETHYL SULFOXIDE; DIMETHYL SULPHOXIDE; DMS-70; DMS-90; DOLIGUR; DOMOSO; DMSO; DROMISOL; GAMAL 90; HYADUR; INFILTRINA; METHANE, SULINYLBIS-; NSC-763; SOMIPRONT; SQ 9453; SYNTEXAN.⁴

ETHANOL


CAS Registry Number: 64-17-5
RTECS Reference Number: KQ6300000
M.W.: 46.08
Chemical Formula: C₂H₆O
Structure:



Synonyms: Absolute Ethanol; Aethanol(GERMAN); Aethylalkohol (GERMAN); Alcohol; Alcohol Anhydrous; Alcohol dehydrated; Alcohol Ethylique(FRENCH); Alcool Etílico(ITALIAN); Algrain; Alkohol(GERMAN); Anhydrol; Cologne Spirit; Cologne spirits (alcohol)(DOT); Etanolo(Italian); Ethanol; Ethanol(DOT); Ethanol 200 proof; Ethyl Alcohol(DUTCH); Ethyl Alcohol Anhydrous; Ethyl Alcohol(DOT); Ethyl Hydrate; Ethyl Hydroxide, Etylowy Alkohol(POLISH); Fermentation Alcohol; Grain Alcohol; Jaysol S; Methylcarbinol; Molasses Alcohol; NA 1170(DOT); NCI-C03134; Potato Alcohol; Spirits of Wine; Spirit; Tecsol.⁴

FORMAMIDE


CAS Registry Number: 75-12-7
RTECS Reference Number: LQ0525000
M.W.: 45.05
Chemical Formula: CH₃NO
Structure:



Synonyms: Carbamaldehyde; Methanamide.⁴

GLYCEROL/FORMAL

CAS Registry Number: 5464-28-8
RTECS Reference Number: JH8390000
M.W.: 104.12
Chemical Formula: C₄H₈O₃
Structure:



GLYCEROL FORMAL

Synonyms: 1,3-dioxan-5-ol and 1,3-dioxolane-4-methanol(mixture);
GF; glycerin-formale; sericosol-n.⁴

HPLC CONDITIONS

DMSO

Column: Microsorb C18 (Rainin, Inst. Co., Woburn, MA)
Mobile Phase: 3% methanol/water
Flow Rate: 2.0 mL/min.
UV: 220nm, .100 absorbance units full scale(aufs)
Injection Volume: 10 μ L

ETHANOL

Column: Dextro-Pak Column (Waters Chromatography
Division, Milford, MA)
Mobile Phase: reagent grade water
Flow Rate: 1.0 mL/min.
Refractive Index: 8X (attenuation)
Injection Volume: 20 μ L

FORMAMIDE

Column: Microsorb C18 (Rainin Inst., Co., Woburn, MA)
Mobile Phase: reagent grade water
Flow Rate: 2.0 mL/min.
UV: 220nm, .100 absorbance units full scale(aufs)
Injection Volume: 10 μ L

GLYCEROL/FORMAL

Column: Dextro-Pak Column (Waters Chromatography
Division, Milford, MA)
Mobile Phase: reagent grade water
flow rate: 1.0 mL/min.
Refractive Index: 8X (Attenuation)
Injection Volume: 20 μ L

PREPARATION OF STOCK AND STANDARD SOLUTIONS

DMSO

DMSO standard solutions were prepared by weighing out separately 15.7, 33.3, 68.5, and 128.1 mg of DMSO in 100-mL volumetric flasks. Reagent grade water was added to the meniscus in each flask to prepare the following standards:

mg DMSO	CONCENTRATION
15.7	157 mg/L
33.3	333 mg/L
68.5	685 mg/L
128.1	1281 mg/L

ETHANOL

Ethanol standard solutions were prepared by weighing out separately 105.5, 215.4, 435.2 and 848.5 mg of ethanol in 100-mL volumetric flasks. Reagent grade water was added to the meniscus in each flask to prepare the following standards:

mg ETHANOL	CONCENTRATION
105.5	1055 mg/L
215.4	2154 mg/L
435.2	4352 mg/L
848.5	8485 mg/L

FORMAMIDE

Formamide standard solutions were prepared by weighing out separately 11.0, 22.3, 43.6 and 123.6 mg of formamide in 100 mL volumetric flasks. Reagent grade water was added to the meniscus in each flask to prepare the following standards:

mg FORMAMIDE	CONCENTRATION
11.0	110 mg/L
22.3	223 mg/L
43.6	436 mg/L
123.6	1236 mg/L

GLYCEROL/FORMAL

Glycerol/Formal standard solutions were prepared by weighing out separately 123.5, 233.2, 423.7 and 847.6 mg of glycerol/formal in 100 mL volumetric flasks. Reagent grade water was added to the meniscus in each flask to prepare the following standards:

mg GLYCEROL/FORMAL	CONCENTRATION
123.5	1235 mg/L
233.2	2332 mg/L
423.7	4237 mg/L
847.6	8476 mg/L

SAMPLE PREPARATION AND HANDLING

Each aqueous sample containing an excess concentration of the appropriate compound was diluted with reagent grade water to obtain a sample concentration within the upper and lower limits of the standard curve. No further sample preparation was required; each sample and standard was injected a minimum of four times.

RESULTS

Figures 1-4 show HPLC chromatograms obtained from aqueous samples of DMSO, ethanol, formamide and glycerol/formal respectively. Each sample peak is symmetrical and well separated from other peaks in each chromatogram.

Calibration curves were constructed by plotting peak areas for all working standards against their concentration. Standard curve data are presented in Table 1. Precision for all four HPLC methods was determined by injecting a low and high concentration sample four times on three separate days. The mean, standard deviation and relative standard deviation are given in Tables 2 and 3.

In order to evaluate method accuracy, recovery studies were conducted on aqueous samples spiked with each compound separately. This was accomplished by taking sample aliquots at two concentration levels and adding a known volume of the appropriate standard. Each aliquot was then analyzed four times to obtain a mean, standard deviation, relative standard deviation and percent recovery (Tables 4 and 5).

DISCUSSION

High performance liquid chromatography (HPLC) was the analytical method of choice for analyzing DMSO, ethanol, formamide and glycerol/formal because it requires almost no sample preparation and the aqueous samples can be injected directly. Since ethanol and glycerol/formal do not exhibit absorbance of ultraviolet (UV) light, a refractive index detector was used in place of the UV detector. The detectors that the HPLC employs (UV and refractive index) can be easily interchanged or connected in series for the appropriate analytical method. The chromatographic run times were optimized to be short in duration because the aqueous samples were clean and free of any interfering peaks. Since high concentrations of solutes were used in the tests (about 1000-7000 mg/L), our lower detection limits could be extended downward by increasing the detector sensitivity if necessary for trace analysis. Calibration curve (DMSO) data for day one shows a slope that is unusually high when compared to days two and three. The calibration curve for day one was replotted to check for any errors in calculations. Both the original and replotted calibration curve data for day one gave the same value for the slope.

CONCLUSION

Relatively rapid and reliable HPLC methods have been developed for the determination of dimethyl sulfoxide, ethanol, formamide and glycerol/formal in aqueous samples. Direct injection of the clean aqueous samples eliminates the need for tedious sample preparation. If needed, the methods will allow for more sensitive detection limits. The calibration curves for all four compounds were linear over the concentration ranges analyzed. Precision and accuracy data for all four HPLC methods are presented in the report.

LITERATURE CITED

1. Ivey, J.P. and P.R. Haddad. 1987. Determination of dimethyl-sulphoxide using ion-exchange chromatography with ultraviolet absorption detection. J. Chromatogr. 391: 309-314.
2. Iwinski, G. and D.R. Jenke. 1987. Determination of alcohols in pharmaceuticals by ion-exclusion chromatography. J. Chromatogr. 392: 397-405.
3. Walbroehl, Y. and J.W. Jorgenson. 1986. Capillary zone electrophoresis of neutral organic molecules by solvophobic association with tetraalkylammonium ion. Anal. Chem. 58:479-481.
4. Registry of Toxic Effects of Chemical Substances. 1985. U.S. Department of Health, Education and Welfare Public Health Service, Center of Disease Control, National Institute for Occupational Safety and Health. Vol. I and II.

DISTRIBUTION LIST

- 4 Commander
U.S. Army Medical Research and Development Command
ATTN: SGRD-RMI-S
Fort Detrick
Frederick, MD 21702-5012

- 2 Commander
U.S. Army Biomedical Research and Development Laboratory
ATTN: SGRD-UBZ-IL
Fort Detrick
Frederick, MD 21702-5010

- 2 Defense Technical Information Center (DTIC)
ATTN: DTIC-FDAC
Cameron Station
Alexandria, VA 22304-6145

- 1 Commander
U.S. Army Biomedical Research and Development Laboratory
ATTN: SGRD-UBZ-P
Fort Detrick
Frederick, MD 21702-5010

- 2 Commandant
Academy of Health Sciences, U.S. Army
ATTN: HSHA-DCD
Fort Sam Houston, TX 78234-6100

TABLE 1. STANDARD CURVE DATA

DMSO			
DAY	SLOPE	R2 VALUE	Y-INTERCEPT
1	25,071	1.0000	4,984
2	16,175	1.0000	4,999
3	16,394	1.0000	4,976

ETHANOL			
DAY	SLOPE	R2 VALUE	Y-INTERCEPT
1	11,854	0.9996	625,627
2	11,988	0.9998	30,312
3	11,634	1.0000	195,572

FORMAMIDE			
DAY	SLOPE	R2 VALUE	Y-INTERCEPT
1	1,816	1.0000	-2,515
2	1,802	1.0000	-2,305
3	1,794	1.0000	-1,373

GLYCEROL/FORMAL			
DAY	SLOPE	R2 VALUE	Y-INTERCEPT
1	18.95	0.9957	-6,560
2	19.25	0.9972	-5,708
3	19.31	0.9986	-6,877

TABLE 2. PRECISION DATA FOR DMSO AND ETHANOL

DMSO

LOW LEVEL

DAY	Mean(mg/L)	S.D.	R.S.D.(%)
1	160	0.43	0.27
2	159	0.24	0.15
3	162	0.29	0.18

HIGH LEVEL

DAY	Mean(mg/L)	S.D.	R.S.D.(%)
1	805	5.34	0.66
2	725	6.64	0.92
3	711	1.37	0.19

ETHANOL

LOW LEVEL

DAY	Mean(mg/L)	S.D.	R.S.D.(%)
1	1,668	105.8	6.34
2	1,475	99.5	6.75
3	1,328	28.7	2.16

HIGH LEVEL

1	6,301	116.3	1.85
2	6,409	233.9	3.65
3	6,299	69.0	1.09

TABLE 3. PRECISION DATA FOR FORMAMIDE AND GLYCEROL/FORMAL

FORMAMIDE

DAY	Mean(mg/L)	LOW LEVEL	R.S.D.(%)
		S.D.	
1	113	0.50	0.44
2	122	0.63	0.52
3	124	0.38	0.31

HIGH LEVEL

1	950	3.01	0.32
2	1003	5.06	0.50
3	991	3.31	0.33

GLYCEROL/FORMAL

Low Level

DAY	Mean(mg/L)	Low Level	R.S.D.(%)
		S.D.	
1	1273	8.25	0.65
2	1213	7.06	0.58
3	1224	13.93	1.14

High Level

1	6209	70.50	1.14
2	6211	23.51	0.38
3	6258	23.12	0.37

TABLE 4. ACCURACY DATA FOR DMSO AND ETHANOL

DMSO

LOW LEVEL

DAY	AMOUNT ADDED (mg/L)	AMOUNT RECOVERED (mg/L)	S.D.	R.S.D. (%)	REC. (%)
1	167	163	1.34	0.82	97.96
2	169	171	2.27	1.33	100.86
3	164	164	0.43	0.26	100.47

HIGH LEVEL

1	922	947	5.28	0.56	102.74
2	943	955	2.77	0.29	101.25
3	931	944	6.73	0.71	101.49

ETHANOL

LOW LEVEL

DAY	AMOUNT ADDED (mg/L)	AMOUNT RECOVERED (mg/L)	S.D.	R.S.D. (%)	REC. (%)
1	1,230	1,326	68.90	5.19	107.83
2	1,230	1,306	36.38	2.79	106.20
3	1,226	1,263	63.22	5.00	103.02

HIGH LEVEL

1	6,148	6,050	50.99	0.84	98.41
2	6,148	6,108	55.21	0.90	99.35
3	6,131	6,118	24.51	0.40	99.79

TABLE 5. ACCURACY DATA FOR FORMAMIDE AND GLYCEROL/FORMAL

FORMAMIDE

LOW LEVEL

DAY	AMOUNT ADDED (mg/L)	AMOUNT RECOVERED (mg/L)	S.D.	R.S.D. (%)	REC. (%)
1	143	142	0.37	0.26	99.29
2	143	140	1.89	1.35	98.18
3	143	141	0.48	0.34	99.11

HIGH LEVEL

1	949	964	4.36	0.45	101.63
2	949	954	4.49	0.47	100.60
3	949	948	2.35	0.25	99.95

GLYCEROL/FORMAL

LOW LEVEL

DAY	AMOUNT ADDED (mg/L)	AMOUNT RECOVERED (mg/L)	S.D.	R.S.D. (%)	REC (%)
1	1331	1266	42.96	3.39	95.09
2	1331	1201	14.59	1.21	90.25
3	1331	1259	12.64	1.00	94.61

HIGH LEVEL

1	6656	6779	20.19	0.30	101.85
2	6656	6635	12.61	0.19	99.69
3	6656	6705	17.91	0.27	100.73